

Alto™ Knee Orthoses

U.S. Orthotics Alto™ Knee Orthoses provide excellent control and support for ligament injuries and instabilities. An ideal brace for post MCL and ACL instabilities as well as medial and lateral instabilities. Our adjustable ROM design features polycentric adjustable hinges for convenient dial-in adjustment of flexion/extension control. Stops at 0°-15°- 30°- 60°- 90° to full ROM. Horizontal latch-top cap cover design ensures there is no inadvertent displacement of hinge pins.



MODEL K-901

13" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.

2 Velcro circumference straps with elastic compression section.



MODEL K-903

17" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.

4 Velcro circumference straps with elastic compression section.



MODEL K-904

Side flexion/extension hinges come fully encased in Neoprene pockets.

13" length with easy pull-up tabs.

2 Velcro circumference straps with elastic compression section.



MODEL K-905

Anterior open wrap-around style for patients who cannot pull on a traditional pull-up support.

13" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.

2 Velcro circumference straps with elastic compression section.

Open popliteal.



MODEL K-907

Anterior open wrap-around style for patients who cannot pull on a traditional pull up support.

17" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.


4 Velcro circumference straps with elastic compression section.

Open popliteal.



ADJUSTABLE HINGE JOINT

Adjustable polycentric hinge joint for flexion/extension control with stops at 0°-15°- 30°- 60°- 90° positions to full ROM.



SIZE CHART

SIZE KNEE CIRCUMFERENCE

Ex-Small.....12"

Small13"

Medium14"

Large.....15"

SIZE KNEE CIRCUMFERENCE

Ex-Large.....16"

2X-Large17"

3X-Large18"

4X-Large.....19"



U.S. Orthotics[®] INC.

800-825-5228

www.usorthotics.com

info@usorthotics.com