



AFO



MODEL B2
Semi Rigid AFO

3/16" natural polypropylene prefabricated semi-rigid Ankle Foot Orthosis. Moderate strength with Medio-Lateral stability. Velcro closure.

Specify left or right.



MODEL B4
Leaf Spring AFO

3/16" natural polypropylene prefabricated leaf spring Ankle Foot Orthosis. Minimal strength. Velcro closure.

Specify left or right.



MODEL B6
AFO with Rolled Calf

3/16" natural polypropylene prefabricated Ankle Foot Orthosis with a rolled calf. Velcro closure.

Specify left or right.



MODEL BSW
Swedish AFO

Leaf spring injection molded polyethylene Ankle Foot Orthosis with open heel and calf. Flexible yet rigid vertical aspect.

SIZE **HEIGHT**
Men's.....12"
Women's.....10 1/2"

MODELS B2 / B4 / B6 SIZING CHART

SIZE	CALF CIRC.	HEIGHT	LONG FOOT
Small.....	10" – 11".....	12 1/2".....	8
Medium.....	11" – 12".....	14".....	8 3/4"
Large.....	11 1/2" – 13 1/2".....	14 1/2".....	9 1/2"
Ex-Large.....	15 1/2" – 17 1/2".....	16".....	10 1/4"

INDICATIONS: Neurological or orthopedic conditions; Lower leg instabilities.



MODEL APS

AFO Strap Kits
Kit consists of: AFO Pad, Reinforced Velcro strap and Chafe tab
Comes complete with an ultrasonically welded strap – with the hook end reinforced with a strong plastic washer, metal link and plastic chafe for easy assembly.

Select Color: Black or White
Select width: 1" – 1 1/2" – 2"



MODEL AP

AFO Pads
Individual plain felt pads – use your own straps & chafe tabs.
Select Color: Black or White
Select pad width: 1" – 1 1/2" – 2"