

ALTO[®] ROM Knee Orthoses



MODEL K-902 L1832

Alto™ ROM Knee Orthoses 13" Pull-On

Alto™ ROM Knee Orthoses Model K-902 is a 13" pull-on hinged knee support. 2 circumference straps and flexion/extension control side hinges. Provides excellent control and support for ligament injuries and instabilities.

Latch-top cap cover design ensures there is no inadvertent displacement of hinge pins. An ideal brace for post MCL and ACL instabilities.

Our adjustable R.O.M. design features polycentric adjustable hinges for convenient dial-in adjustment of flexion/extension control.

ROM Stops at 0° – 15° – 30° – 60° – 90° to Full ROM

NEOPRENE or DRY FIT Fabric

Sizes	Knee Circumference	Indications
Ex-Small.....	12"	Post-operative management & rehabilitation. Locked or limited motion control of knee. Mild to moderate ACL, PCL, MCL, LCL. Combined ligament instabilities. Medial and lateral collateral ligament instabilities.
Small.....	13"	
Medium	14"	
Large	15"	
Ex-Large	16"	
2X-Large.....	17"	
3X-Large.....	18"	
4X-Large.....	19"	

Customs upon request.

