

# Shoulder Immobilizers



**MODEL SM-150-U**  
**DELUXE SHOULDER IMMOBILIZER**

Soft, perforated foam sling and detachable swathe. Velcro® closures. Sling fits over each shoulder (not around neck) with Velcro® front height adjustment straps to position patient's arm. The 4 1/2" wide swathe is adjustable in circumference to keep arm securely against the body. Perforated for coolness.

Universal.  
One size fits most adults.  
52" maximum circumference.  
Fits both left and right.



**MODEL SM-155-U A4565**  
**SHOULDER IMMOBILIZER**

Body band with an adjustable 4" upper arm cuff and over-the-shoulder adjustable wrist and hand cuff. All foam lined material with Velcro® closures. Perforated for coolness.

Universal.  
One size fits most adults.  
46" maximum circumference.  
Fits both left and right.



**MODEL SM-156-U A4565**  
**SHOULDER IMMOBILIZER**

Body band with a 4" adjustable upper arm cuff and 3" adjustable wrist cuff. Arm and wrist cuffs are held in position with Velcro®. Perforated for coolness.

Universal.  
One size fits most adults.  
46" maximum circumference.  
Fits both left and right.