

# Spina<sup>®</sup> II Spinal Orthosis Model 1200

CLASSIC LSO



## SPINA<sup>®</sup> II LIGHTWEIGHT PLASTIC SPINAL ORTHOSIS SOFT FRONT L0633

Traditional plastic chairback orthoses lined with plush soft-foam.  
Lightweight & durable. The 2" straps are made of double-sided Velcro  
webbing, adjustable at both ends which provide a wide range of flexibility.

12" back height and 9" soft front.

Indications: Chronic Back Pain, Compression Fracture, Herniated / Bulging  
Disk, Laminectomy, Ligament Strain, Muscle Sprain, Spinal Stabilization,  
Spinal Stenosis, Spondylolisthesis, Spondylolysis

SIZE	HIP CIRCUMFERENCE
Small.....	30" – 35"
Medium.....	36" – 41"
Medium.....	36" – 41"
Large.....	42" – 47"
Ex-Large.....	48" – 52"
2X-Large.....	52" – 56"