

Spina[®] III Spinal Orthosis

CLASSIC LSO



MODEL 1219 PDAC APPROVED L0627

Rigid anterior panel and soft-foam lined posterior rigid frame

Traditional chair back orthosis with rigid anterior panel and soft-foam lined posterior rigid frame. The 2" double-sided Velcro closure straps are adjustable at both ends which provide a wide range of flexibility. Easy cross-pull straps provide even compression.

12" back height

Standard 9" front with a 7" rigid anterior insert. Model 1219

Low profile 7" front with a 5" rigid anterior insert. Model 1218

SIZE	HIP CIRCUMFERENCE
Small.....	30" – 35"
Medium.....	36" – 41"
Large.....	42" – 47"
Ex-Large.....	48" – 51"
2X-Large.....	52" – 55"