

# Ostomy Belts




## MODEL OSB

### Standard Hole Placement

Elastic ostomy belts with reinforced stoma opening. Thin enclosed polypropylene panel provides additional support to the abdomen around the stoma. Velcro closure.

Short or long  
Black or white  
Belt Height: 4" – 6" – 8" – 10" – 12"  
Access Hole Size: 2.5" – 3" – 3.5" – 4"

### Standard access hole location – centered on belt height



SIZE	CIRCUMFERENCE IN LINE WITH STOMA
Short	28" – 38"
Long	38" – 48"


## MODEL OSB-C

### Custom Hole Placement

Elastic ostomy belts with reinforced stoma opening. Thin enclosed polypropylene panel provides additional support to the abdomen around the stoma. Velcro closure.

Short or long  
Black or white  
Belt Height: 4" – 6" – 8" – 10" – 12"  
Access Hole Size: 2.5" – 3" – 3.5" – 4"

### Custom access hole location



SIZE	CIRCUMFERENCE IN LINE WITH STOMA
Short	28" – 38"
Long	38" – 48"
Custom	_____

### CUSTOM OSB-C STEPS

1. Select belt length: short / long / custom
2. Select patient Left or Right
3. Select belt color: Black or White
4. Select belt height: 4" – 6" – 8" – 10" – 12"
5. Select access hole size: 2.5" – 3" – 3.5" – 4"
6. Positioning the custom access hole:
  - A. From the edge of the belt, to the bottom of the access hole.
  - B. From the edge of the belt end, to the side of the access hole.

Fillable form can be found on OSB page at [www.usorthotics.com](http://www.usorthotics.com)