

Z-Cross™ Lumbosacral Orthoses



MODEL Z4-U L0627 / L0642



MODEL Z5-U L0631 / L0648



MODEL Z8-U L0637 / L0650

U.S. Orthotics Z-Cross™ Spinal Orthoses are made with a 6:1 mechanical advantage compression strapping system, utilizing dual pull string pulley straps. One smooth pull securely forms the Z-Cross™ to the patient's torso and provides intra-abdominal compression, reducing the load on inter-vertebral discs.

- * AP Control LSO
- * Bilateral pull handles
- * 2 height posterior panels 12" & 14"
- * Conforms to any body shape
- * Smooth dual-pull pulley compression system
- * Universal Sizing - One size fits most.
- * Fits 28" - 48" with extensions up to 62"

PDAC Approved
Model - Custom Fit / Off the Shelf (OTS)

Z4 L0627 / L0642
Z5 L0631 / L0648
Z8 L0637 / L0650



Anterior View



Lateral Panels



Extension Panels