

Soft Cervical Collars



MODEL CU L0120

UNIVERSAL CIRCUMFERENCE

1" soft poly foam cervical collar. Cotton stockinette cover. Velcro® closure.

HEIGHTS 2" - 2.5" - 3" - 3.5" - 4"

Universal circumference fits approximately 15" - 21".



MODEL CC L0120

SIZED BY CIRCUMFERENCE

1" soft poly foam cervical collar. Cotton stockinette cover. Velcro® closure.

HEIGHTS 2" - 2.5" - 3" - 3.5" - 4"

Small..... 12" X 14"

Medium..... 15" X 17"

Large..... 17" X 19"

Ex-Large..... 19" X 21"