

FULLBACK[®] XPA Lumbosacral Support



MODEL FULLBACK[®] XPA PDAC APPROVED L0627

Back Support with Anterior & Posterior Inserts

Single piece designed elastic back support with anterior and posterior thermo-formable plastic inserts. Top tension strap for additional pressure. Latex-free elastic. Pre-molded Kydex padded inserts included.

Black or White

3 Models

- 8" back height tapered to a 4" front.
- 10" back height tapered to a 6" front.
- 12" back height tapered to a 6" front.

MODEL XPA

- Small..... 29" – 33"
- Medium..... 33" – 37"
- Large 37" – 41"
- Ex-Large..... 41" – 45"
- 2X-Large..... 45" – 49"
- 3X-Large..... 49" – 53"
- 4X-Large..... 53" – 57"
- 5X-Large..... 57" – 61"