

Anterior Opening Neoprene Knee Support



MODEL KW6 PDAC APPROVED L1810

Hinged wrap-around knee support

12" wrap-around neoprene knee support with polycentric side hinges and 2" adjustable Velcro straps. Excellent for patients who cannot pull on a traditional neoprene knee support.

Open popliteal. Opens flat. Wraps around. Easy to apply.

SIZE	KNEE CIRCUMFERENCE
Small.....	11" – 13"
Medium.....	13" – 15"
Large	15" – 17"
Ex-Large	17" – 19"
2X-Large.....	19" – 21"
3X-Large.....	21" – 23"