

Alto[®] Knee Orthoses

ALTO[®]
KNEE ORTHOSIS

Alto[®] Knee Orthoses line provides excellent control and support for ligament injuries and instabilities. An ideal brace for post MCL and ACL instabilities as well as medial and lateral instabilities. Our adjustable ROM design features adjustable hinges for convenient dial-in adjustment of flexion/extension control. Latch-top cover design ensures there is no inadvertent displacement of hinge pins. Available in dry-fit fabric or standard Neoprene.



MODEL K-901
With Condyle Pads

13" pull-on hinged knee support



MODEL K-902
Without Condyle Pads

13" pull-on hinged knee support



MODEL K-903
17" Pull-On

17" pull-on hinged knee support



MODEL K-904
Covered Hinge

13" pull-on hinged knee support



MODEL K-905
Wrap Around

13" anterior opening wrap-around hinged knee support with condyle pads



MODEL K-906
Wrap Around

13" anterior opening wrap-around hinged knee support without condyle pads



MODEL K-907
17" Wrap Around

17" anterior opening wrap-around hinged knee support



ADJUSTABLE HINGE JOINT

Adjustable hinge joint for flexion/extension control with stops at:
 0°-15°- 30°- 60°- 90° positions
 to full ROM

WRAP-AROUND SIZE CHART

KNEE CENTER

Small	11"-13"	2XL	19"-21"
Medium	13"-15"	3XL	21"-23"
Large	15"-17"	4XL	23"-25"
Ex-Lg	17"-19"		

PULL ON SIZE CHART

KNEE CENTER

Ex-Small	12"	Ex-Lg	16"
Small	13"	2XL	17"
Medium	14"	3XL	18"
Large	15"	4XL	19"